

Identity Testing of Raw Materials using Raman Spectroscopy: Active Pharmaceutical Ingredients

Relevant for: Pharma

Cora 5001 together with the pharma compliant Spectroscopy Suite desktop software is a fast tool for identity testing of active pharmaceutical ingredients (APIs). Method development and validation according to USP regulation is easily possible.



The pharmaceutical production industry is one of the most heavily regulated markets in the world. From the delivery of raw materials to the packaging of final products, the production processes need to be monitored according to legal requirements ensuring product quality and minimizing risks for the patients in all parts of the world. The Active Pharmaceutical Ingredients (APIs), the pharmacologically active part of a drug product, need to be tested among other things for their chemical identity as stated by international pharmacopeias and GMP guidelines.[1], [2], [3] Here, fast and accurate methods are required to quickly assess whether a batch is fit for further use or needs to be rejected. Anton Paar's Cora 5001 Raman spectrometers together with the Spectroscopy Suite desktop software offer fast and reliable ID testing using algorithmic database matching. Method development and validation are straightforward and processes are easy to handle due to the intuitive user interface. As the software is fully compliant to CFR 21 part 11 according to FDA, data integrity and a complete audit trail are ensured.[4]

1 Active Pharmaceutical Ingredients (APIs)

According to GMP an active pharmaceutical ingredient (API) or drug substance is "any substance or mixture of substances intended to be used in the manufacture of a drug (medicinal) product and that, when used in the production of a drug, becomes an active ingredient of the drug product".[4] Furthermore, GMP states that

these "substances are intended to furnish pharmacological activity or other direct effect in the diagnosis, cure, mitigation, treatment, or prevention of disease or to affect the structure and function of the body."[4]

According to international pharmacopeias identity testing is required for every API used in drug production as stated in the respective chapters.[5], [6]

2 Requirements for ID Testing with Raman Spectroscopy according to USP

According to US Pharmacopeia chapter 197, Raman spectroscopy may be used for identity testing of raw materials such as APIs.[6] The following recommendations apply: measurement device should be prepared for the measurement according to the manufacturer's recommendations and a performance check should be conducted prior to use. Spectra of sample and references should be collected using the same instrumental parameters. The resulting spectra are then compared qualitatively. Where available a USP reference standard should be used as reference sample.[5] The qualitative comparison can be conducted algorithmically using chemometric models.[7] The method has to be validated according to USP 1225.[8] For qualitative analyses the specificity and the robustness of the method have to be demonstrated during method validation.[7], [8]



Figure 1: Anton Paar's Cora 5001 Raman spectrometer series (left: Direct version; right: Fiber version).

Table 1: List of measured APIs.

API	Drug group	abbreviation
Acetaminophen	Analgesic, antipyretic	AAP
Amoxicillin	antibiotic	AMX
Aspirin	Analgesic, anti-inflammatory drug, antipyretic	ASP
Azithromycin	antibiotic	AZM
Caffeine	CNS stimulant	CAF
Ciprofloxacin	antibiotic	CFX
Clarithromycin	antibiotic	CLM
Dopamine HCl	Anti-hypotensive, treatment of bradycardia	DOP
Hyaluronic acid sodium salt	glycosaminoglycan	HYA
Hydrocortisone	Corticosteroid, Immunosuppressive drug	HYD
Ibuprofen	Non-steroidal anti-inflammatory drug, anti-pyretic, analgesic	IBU

3 ID Testing using Raman Spectroscopy

3.1 Experimental Setup

For ID Testing of the APIs a Cora 5001 Direct Raman spectrometer with an excitation wavelength of 785 nm was used (see Figure 1). The instrument is compliant to USP 858 and therefore perfectly suited for ID testing using Raman spectroscopy.[9] It was controlled via Anton Paar's Spectroscopy Suite desktop software, which fully complies to CFR 21 part 11 according to FDA, ensuring data integrity.[4]

3.2 Method Development

The instrument was checked for correct calibration with the supplied System Suitability Test (SST) Standard prior to use. As a first step to method development, test runs for optimizations of measurement conditions were conducted. The implemented autofocus and automatically determined exposure time render this process swift and hassle-free. The obtained parameters were saved in a method which was submitted, reviewed and approved by the responsible users. Subsequently, a good quality library spectrum of each reference substance was measured using the previously developed method and submitted for approval. After approval by the responsible lab manager the spectra were fit to be used in a spectral library.

3.2.1 Spectral Library Development, Publication and Approval

For each substance a library entry with the corresponding spectrum was created. The software

offers the possibility to enter additional information on the substance such as further names or GHS symbols as well as a comment field for other relevant information (e.g. CAS number, purity, supplier). Subsequently, each substance entry and the whole library are submitted for approval. After approval the library is fit for use in an ID testing method.

3.2.2 ID Testing

An identification method is created using the same parameters as used for the measurements of the library spectra, as required by USP. The software offers to copy the method used to acquire the reference data and to add identification parameters for the intended analysis, therefore equality of the measurement parameters can be easily achieved. The new method is signed (submit, review, approve) as well by the responsible employees. Subsequently, samples of each API (see Table 1) were measured and identified.

3.3 HQI Matching Algorithm

For the comparison of reference and sample spectra the Anton Paar Spectroscopy Suite offers an HQI matching algorithm. The results of the algorithmic matching lie between zero for no correspondence and 100 for identical spectra.

By defining a suitable HQI threshold the identification or verification of substances can easily be achieved.

4 Results and Method Validation

Table 2 shows the resulting matches of each sample against each reference spectrum using the HQI (Hit Quality Index) matching algorithm. It is clearly visible that the substances are identified correctly with an HQI threshold of 94. The specificity of the method is clearly visible when comparing the spectra of structurally very similar substances such as Azithromycin and Clarithromycin. The HQI difference is > 12 and therefore the substances are still well separable with the implemented algorithm though they are chemically closely related.

Other Antibiotics which are not structurally similar are well separated with HQI differences above 87. The visual comparison of the different spectra of the antibiotics as shown in Figure 2 again illustrates the discrimination power of Raman spectroscopy.

Even clearer spectra are obtained for typical APIs of over-the-counter (OTC) analgesic drugs such as ibuprofen, acetyl salicylic acid, or acetaminophen (Figure 3). The spectra are less crowded due to the fact that the molecules which were compared are relatively small. The HQI results differ by more than 81 again showing an excellent discrimination using the simple HQI algorithm.

Table 2: HQI Matrix showing specificity and robustness (average of triple measurements) of ID testing method. Marked in red are antibiotics, marked in yellow are over-the-counter drugs.

		Test samples (average of 3 measurements)										
		AAP	ASP	AMX	AZM	CAF	CFX	CLM	DOP	HYA	HYD	IBU
Reference substances	AAP	99,98	11,88	12,86	2,52	4,83	13,53	1,95	5,99	2,21	35,60	8,42
	ASP	15,17	97,01	2,62	0,00	2,53	12,54	0,04	13,44	0,02	44,10	11,89
	AMX	13,03	2,05	99,98	0,13	0,01	3,59	0,44	1,47	0,18	0,55	8,17
	AZM	3,09	0,01	0,05	98,31	2,19	11,01	69,88	4,10	60,23	0,21	7,48
	CAF	4,99	3,32	0,01	1,99	99,88	0,25	1,32	5,83	1,23	1,78	0,88
	CFX	13,44	10,51	3,56	9,25	0,21	99,99	5,47	4,75	11,93	9,45	2,15
	CLM	2,22	0,06	0,28	85,28	1,45	7,69	96,53	4,37	49,00	0,10	15,36
	DOP	5,43	14,56	1,41	2,75	5,59	4,19	3,50	99,78	1,37	3,37	8,12
	HYA	1,84	0,10	0,20	48,27	0,85	11,08	34,01	1,92	99,35	0,19	2,45
	HYD	39,04	37,32	0,56	0,90	2,44	10,95	0,50	4,98	1,13	97,84	7,37
	IBU	7,90	11,33	8,17	7,68	0,77	2,03	17,72	8,47	2,46	5,86	99,95




Figure 2: Raman spectra of four different antibiotics as shown in the processing menu of the Anton Paar Spectroscopy Suite software. Spectra have been normalized to the highest peak.




Figure 3: Raman spectra of typical APIs in OTC analgesics as shown in the processing menu of the Anton Paar Spectroscopy Suite software. Spectra have been normalized to the highest peak.

Measurements were performed in triplicate to also evaluate repeatability, standard deviations varied between 0,00004 to 1,36 of the HQI. It is recommended to validate a method by measuring spectra on different days and by different instrument operators.

After validation the method is fit for use in the production.

5 Conclusion

Raman spectroscopy using Anton Paar's Cora 5001 Raman spectrometer together with the Spectroscopy Suite Desktop Software speeds up identity testing for APIs. Validation of the method demonstrates specificity and robustness of the procedure and is established equally fast as the test measurements themselves.

6 References

- [1] USP. Oral Drug Products – Product Quality Tests <2> In: USP-NF. Rockville, MD: USP; Oct 12, 2021.
DOI: https://doi.org/10.31003/USPNF_M3211_05_01
- [2] USP. Topical and Transdermal Drug Products – Product Quality Tests <3>. In: USP-NF. Rockville, MD: USP; Oct 12, 2021.
DOI: https://doi.org/10.31003/USPNF_M4033_05_01
- [3] 21 CFR part 210. Current Good Manufacturing Practice in Manufacturing, Processing, Packing, or Holding of Drugs; General. U. S. Food and Drug Administration, retrieved Nov 11, 2021.
<https://www.ecfr.gov/current/title-21/chapter-I/subchapter-C/part-210>
- [4] 21 CFR part 11. Electronic records; electronic signatures. U. S. Food and Drug Administration, retrieved Oct 12, 2021.
<https://www.accessdata.fda.gov/SCRIPTS/cdrh/cfdocs/cf/cfr/CFRSearch.cfm?CFRPart=11&showFR=1>
- [5] Ec.europa.eu. 2015. EudraLex - Volume 4 Good manufacturing practice (GMP) Guidelines. - European Commission. [online] Available at: https://ec.europa.eu/health/medicinal-products/eudralex/eudralex-volume-4_en, Retrieved Nov 8, 2021.

- [6] USP. Identification Tests – General <191>. In: USP-NF. Rockville, MD: USP; Oct 12, 2021.
DOI: https://doi.org/10.31003/USPNF_M98940_09_01
- [7] USP. Spectroscopic Identification Tests <197>. In: USP-NF. Rockville, MD: USP; Oct 12, 2021.
DOI: https://doi.org/10.31003/USPNF_M98947_05_01
- [8] USP. Chemometrics <1039>. In: USP-NF. Rockville, MD: USP; Oct 26, 2021.
DOI: https://doi.org/10.31003/USPNF_M2345_02_01
- [9] USP. Validation of Compendial Procedures <1225>. In: USP-NF. Rockville, MD: USP; Oct 26, 2021.
DOI: https://doi.org/10.31003/USPNF_M99945_04_01
- [10] USP. Raman Spectroscopy <858>. In: USP-NF. Rockville, MD: USP; Oct 12, 2021.
DOI: https://doi.org/10.31003/USPNF_M8188_02_01

Contact Anton Paar GmbH

Tel: +49 511 40095-0 | www.anton-paar.com
application-optotec@anton-paar.com